



Self-Retaining Soft Tissue Retractors

Used for soft tissue retraction during hip & shoulder surgery



PRODUCT NO'S:	
1571	Overall Length: 7.75"
Included Modular Snap-In Blades	
1570-01 [Double Blunt-tipped]	Blade Depth: .75"
1570-02 [Rounded]	Blade Depth: 1"
1570-03 [Chandler-type]	Blade Depth: 1.75"

Center Retraction Blade No Longer Available



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

INNOMED

1.800.548.2362