

# APC Proximal Femoral Elevator

*Elevates the proximal femur during total hip or hemi-arthroplasty surgery*



Designed to elevate the proximal femur during total hip or hemi-arthroplasty surgery. Its unique design provides excellent access to the intramedullary canal. The elevator's geometry incorporates serrated edges to grip and elevate the proximal femur.

Designed by  
APC, Inc.

Product No's  
3421-00 [Standard]  
Overall Length: 10.75"  
Blade Width at Widest: 63mm

3421-01 [Small]  
Overall Length: 10.75"  
Blade Width at Widest: 50mm

**FREE TRIAL**  
ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION



©2009 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**  
912 236 0000  
Fax 912 236 7766  
Internet [www.innomed.net](http://www.innomed.net)  
eMail [info@innomed.net](mailto:info@innomed.net)

**Innomed**  
ORTHOPEDIC INSTRUMENTS  
103 Estus Drive Savannah GA 31404

