

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2000 Innomed, Inc., All Rights Reserved.

Innomed
Instrument Division
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Acetabular Graft Protectors

The anatomically shaped and sized acetabular graft protector allows suction of blood after acetabular bone grafting. The protector helps protect the grafted bone from being dislodged during suctioning. Available in two sizes.

Designed by
Steven B. Zelicof, MD

Product No's: 1690-44 [44mm]
1690-52 [52mm]

*Helps protect grafted bone from
being dislodged during suctioning*

