


# Acromioplasty Retractor

*Designed to retract and protect the humeral head during resection of the inferior acromial surface*

The two prongs hook the posterior aspect of the acromion for retraction while the file is used to smooth rough edges of the acromion post-resection

Product No: S3008  
Overall Length: 9"  
Blade Width: 18mm



**FREE TRIAL**  
ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

**Phone** 800 548 2362  
912 236 0000  
**Fax** 912 236 7766  
**Internet** [www.innomed.net](http://www.innomed.net)  
**eMail** [info@innomed.net](mailto:info@innomed.net)

**Innomed**  
ORTHOPEDIC INSTRUMENTS  
103 Estus Drive Savannah GA 31404

