



Amstutz Hip Resurfacing Instruments

Designed by Harlan C. Amstutz, MD

*Amstutz Femoral Head-Neck Elevator
Used for acetabular exposure*



PRODUCT NO'S:
3410 [Wide Blade] Overall Length: 12.25" Blade Width at Widest: 67mm
3410-01 [Narrow Blade] Overall Length: 12.25" Blade Width at Widest: 50mm



*Amstutz Femoral Retractor
Used for acetabular exposure*

PRODUCT NO:
6410 Overall Length: 10.75" Blade Width: 18mm



*Amstutz Inferior
Acetabular Retractors*

PRODUCT NO'S:
7690-01 [Right] Overall Length: 8.375" Blade Width: 54mm
7690-02 [Left] Overall Length: 8.375" Blade Width: 54mm

Cement Packer & Trimmer

PRODUCT NO:
4995 Overall Length: 9.75"



*Amstutz Acetabular
Exposure Pin System*

PRODUCT NO'S:	
1200-00 [Set: Ins/Ext & Two Pins]	
1200-0A [Set: Ins/Ext & Two Pins w/Stop]	
Also sold Individually:	
1200-01 [Inserter/Extractor]	
1200-02 [Pin] Overall Length: 4.5" Pin Depth: 2" Pin Diam.: 3.8mm	1200-03 [Pin with Stop] Overall Length: 4.5" Pin Tip-to-Stop Depth: .75" Pin Diameter: 3.2mm

*Amstutz Angled Mirrors
Used for viewing hard-to-see areas
around the acetabulum*



PRODUCT NO'S:	
4860 [Left] Overall Length: 8" Mirror #1: 32x13mm Mirror #2: 19x16mm	4861 [Right] Overall Length: 8" Mirror #1: 32x13mm Mirror #2: 19x16mm

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74



1.800.548.2362