

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2001 Innomed, Inc., All Rights Reserved.

Innomed
Instrument Division
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Angled Lowman-Type Bone Clamp

The offset distance between the jaws and handle of the clamp allow space for free and easy access to use a drill or screwdriver. The angled clamp and more-open and thinner jaws facilitate easier use in deep incisions. The angled shaft also acts as a self-retaining retractor. The tightening handle is scalloped to lessen slippage when tightening or untightening.

Designed by
John J. McLeod, Jr., M.D.

Product No: 1770



*Angled for easier insertion of
the jaws around the bone*

