



Bowen Arthroscopic Shoulder Elevator

Used to enlarge the Bankart lesion

Permits access to the scapular neck. Terminal bend matches the angle of the scapular neck from a standard anterior portal.

Product No: 2305
Overall Length: 9.625"
Handle Length: 3.5"



Bowen Arthroscopic Shoulder Rasp

Used to abrade the anterior scapular neck to stimulate a vascular healing response

Terminal bend matches the angle of the scapular neck from a standard anterior portal.

Product No: 2310
Overall Length: 9.625"
Handle Length: 3.5"

FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

