

# Baldwin Lateral Soft Tissue Retractor

*Designed to hold back the fat pad and soft tissues during total knee arthroplasty*

The fenestrated paddle helps hold back the fat pad and soft tissues, while the two sharp-tipped prongs help penetrate the soft tissue, but have flat surfaces that rest against the side of the tibia and help prevent rotation of the instrument.

Designed by James L. Baldwin, MD

**PRODUCT NO:**

6312

Overall Length: 9.875"

Pad Dimensions: 38mm x 15mm

Prong Depth: 22mm

MADE  
IN USA



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

© 2012 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

**INNOMED**

1.800.548.2362