


Barger Anterior Capsular Retractor

Designed by Daniel Barger, PA-C and Ryan Pokorney, DO

Designed for retraction of the medial capsule in anterior total hip arthroplasty



PRODUCT NO:

6423

Overall Length: 12" (30,5 cm)

Blade Depth: 3.75" (9,5 cm)

Blade Width: 1.3" (3,3 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362