

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2003 Innomed, Inc., All Rights Reserved.

Innomed
 Instrument Division
 103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS
 CALL FOR MORE INFORMATION

Bargo Bone Holding Clamp

Designed to aid in the reduction of various fractures such as: spiral, transverse, compound, oblique, or butterfly. The clamp can also be used to secure a plate in place while the screw holes are being drilled and screws inserted. The fracture site can also be manipulated with the clamp being used as a lever. Available in two sizes, large and small, it has teeth in the jaws for a better grip and a ratchet locking handle for use on various bone diameters.

Designed by
 Lonnie Bargo, CST/CFA

Product No's:
 1895-02 [Large]
 1895-01 [Small]

*Designed to aid in the reduction of various fractures,
 and can help secure a plate in place during installation*

