

# Bechtold Ergonomic Orthopedic Mallet

Designed by Dustin Bechtold, MD

*Ergonomically designed for forward and backward strikes, featuring an ergonomic handle with a tamp*



- ▶ Stainless steel head and shaft with an aluminum handle with a right-handed grip
- ▶ Large and small striking heads with smooth surface
- ▶ Palmar side of the mallet features a flat surface to slide along a broach or impacting type instrument for back slapping and serves well as an additional striking surface

#### PRODUCT NO:

7822  
Weight: 2.7 lbs (1.22 kg)  
Overall Length: 10.75" (27.3 cm)  
Head Width: 4" (10.2 cm)  
Large Head Diameter: 2" (5.1 cm)  
Small Head Diameter: 1.5" (3.8 cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.