


# Bent Hohmann Retractors—Narrow

*Helps retract tissues at the margins of the joint*

Useful for retracting tissues at the margins of the joint. Can be passed over the margins of the joint and held in place with weights or by hand.

The OrthoLucent™ version is made of a lightweight carbon fiber material, which is completely radiolucent, helps to prevent from scratching component surfaces, and can be steam sterilized.



PRODUCT NO'S:
<b>7110 [Narrow]</b> Overall Length: 9.75" Handle Length: 7" Blade Width: 19mm Depth from Bend: 4.75"
<b>7110-R [OrthoLucent™ Narrow]</b> Overall Length: 9.75" Handle Length: 7" Blade Width: 19mm Depth from Bend: 4.75"
<b>7110-01 [Narrow with Extra Long Handle]</b> Overall Length: 11.5" Handle Length: 10" Blade Width: 19mm Depth from Bend: 4.75"
<b>7115 [Short-tipped Narrow]</b> Designed by Carl DiRaimondo, MD Overall Length: 9.25" Handle Length: 7" Blade Width: 19mm Depth from Bend: 4.25"
<b>7115-01 [Short-tipped Narrow with Extra Long Handle]</b> Overall Length: 11" Handle Length: 10" Blade Width: 19mm Depth from Bend: 4.25"

OrthoLucent™ is a trademark of MedX Composites

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
 www.innomed.net info@innomed.net

912.236.0000 Phone  
 912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73  
 Fax +41 41 740 67 74



1.800.548.2362