

# Bent Hohmann Retractors—Narrow

*Helps retract tissues at the margins of the joint*

Useful for retracting tissues at the margins of the joint. Can be passed over the margins of the joint and held in place with weights or by hand.

*The OrthoLucent™ version is made of a strong, lightweight carbon fiber PEEK composite material, which is completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.*

## PRODUCT NO'S:

7110 [Narrow]
Overall Length: 9.75" (23,8cm)
Handle Length: 7" (17,8cm)
Blade Width: 19mm
Depth from Bend: 4.75" (12,1cm)
7110-R* [OrthoLucent™ Narrow]
Overall Length: 9.75" (23,8cm)
Handle Length: 7" (17,8cm)
Blade Width: 19mm
Depth from Bend: 4.75" (12,1cm)
7110-01 [Narrow with Extra Long Handle]
Overall Length: 11.5" (29,2cm)
Handle Length: 10" (25,4cm)
Blade Width: 19mm
Depth from Bend: 4.75" (12,1cm)
7111 [Narrow with Extra Grip Tip]
Overall Length: 9.75"
Blade Width: 19 mm
Depth from Bend: 4.25"
7115 [Short-tipped Narrow]
Overall Length: 8.625" (21,9cm)
Handle Length: 7" (17,8cm)
Blade Width: 19mm
Depth from Bend: 4.4" (11,2cm)
7115-01 [Short-tipped Narrow with Extra Long Handle]
Overall Length: 11" (27,9cm)
Handle Length: 10" (25,4cm)
Blade Width: 19mm
Depth from Bend: 4.25" (10,8cm)
7115-03 [Extra Deep Short-tipped Narrow]
Overall Length: 12.125" (31,1cm)
Handle Length: 9.75" (24,8cm)
Depth from Bend: 6.25" (15,9cm)
Blade Width: 19mm




Short-tipped designed by Carl DiRaimondo, MD

Extra Grip Tip design modification by Alfred A. Durham, MD



## Extra Grip Tip

*Helps to stabilize the retractor around the bone*



Narrow  
#7110

OrthoLucent™ Narrow  
#7110-R

Narrow with Extra Long Handle  
#7110-01

Narrow with Extra Grip Tip  
#7111

Short-tipped Narrow  
#7115

Short-tipped Narrow with Extra Long Handle  
#7115-01

Extra Deep Short-tipped Narrow  
#7115-03

[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

# INNOMED