



Beyer Wedges

Can be used during total knee surgery as a retractor or for assessing intraoperative soft tissue balance, or can be used in same-size pairs to assess flexion/extension balance

Made of delrin to prevent from scratching implant surfaces.

Wider sizes also useful in revision knee surgery.

PRODUCT NO'S:

1204-08 [8mm]
Overall Length: 5.5"

1204-12 [12mm]
Overall Length: 5.5"

1204-16 [16mm]
Overall Length: 5.5"

Designed by Craig A. Beyer, MD



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362