

Modified Blunt Hohmann Retractor

Used for soft tissue retraction

The OrthoLucent™ version is made of a lightweight carbon fiber material, which is completely radiolucent, helps to prevent from scratching component surfaces, and can be steam sterilized.



PRODUCT NO'S:
4550 [Standard] Overall Length: 10.75" Blade Width: 11mm
4550-R [OrthoLucent™] Overall Length: 10.75" Blade Width: 11mm

OrthoLucent™ is a trademark of MedX Composites

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

INNOMED

1.800.548.2362