

# Bolanos Hip Retractor

Designed by Alberto Bolanos, MD


*Designed to help improve femoral exposure while reducing the risk of trochanteric fractures in both posterior and anterior approaches*

Designed to lever against cortical bone at the subtrochanteric level, which helps spare the greater trochanter from stresses. The wrap-around design helps the assistant to retract the entire gluteus medius during femoral preparation.

#### PRODUCT NO'S:

6350-R [Right]  
Overall Length: 12.5"  
Handle Length: 4.5"

6350-L [Left]  
Overall Length: 12.5"  
Handle Length: 4.5"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362