

# Bolanos Shoulder Retractor

Designed by Alberto Bolanos, MD

*Designed for mini-open rotator cuff repairs and shoulder arthroplasty, the contour matches the humeral head and the rounded edge helps avoid trauma to surrounding musculature*

Depth matches girth of most patients, while the comfortable handle makes it easier for assistants to hold.



**PRODUCT NO:**

3222

Overall Length: 7.5" (19,1 cm)

Blade Width at Widest: 1" (2,54 cm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

[www.innomed.net](http://www.innomed.net)



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

# INNOMED

**1.800.548.2362**