

Depth Gauge

Designed to help provide measurement of the depth/length of any bone hole for proper screw length determination

PRODUCT NO:

8015

Overall Length: 7.625" (19,4 cm)
Scale: From 0 to 48 mm



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362