



Curved Cobra Retractors for Total Hip Resurfacing

Left and right retractors can be used with posterior or lateral approach for total hip resurfacing

Placed over the anterior acetabular wall for supero-anterior translation of the femoral head and neck during acetabular preparation. The two prongs pierce the capsule over the anterior wall, then the curve and "twist" of the retractor allows for gentle retraction on the femoral neck.

Designed by
Henry Boucher, MD

Product No's:
6110-01 [Right]
6110-02 [Left]



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

