



# Browner Wire Tightener

Designed by Bruce D. Browner, MD

*Wire is passed through the distal arm hole and into the separate drum holes, and can then be tightened and rotated before being cut with a wire cutter*



**PRODUCT NO:**  
8251  
Overall Length: 6" (15,2 cm)  
Width: 3.75" (9,5 cm)  
Wire Hole Diameters: .125" (3,2 mm)

**PROUDLY  
MADE  
IN THE  
USA**

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**

# INNOMED