

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Calcaneal Spreader

This device has tongs with large surface area, which easily separates the calcaneal osteotomized bone for placement of tricortical bone graft. The large surface area helps prevent the compression of soft calcaneal cancellous bone that is often seen with narrower devices.

Designed by
Michael Forness, DO

Product No's:
1880 [Smooth Pads]
1881 [Grooved Pads]

*Separates the calcaneal osteotomized bone
for placement of tricortical bone graft*

