


# Cannestra Extended Femoral Osteotomy Guide

Designed by Vince Cannestra, MD

Designed for femoral shaft osteotomies during hip revision surgery



Set includes a sterilization case, and a package of 10 short and 10 long threaded bone pins with collar.

#### PRODUCT NO'S:

1158 [Set with Guide, 10 of Each Pin Size]

#### Individual/Replacement Parts:

1158-01 [Osteotomy Guide Only]

Overall Length: 8"

Cutting Slot Length: 7.5"

1225-L [Long Threaded Bone Pins] **Pkg of 10**

Overall Length: 40mm

Collar-to-Tip Length: 25mm

Diameter: 3.2mm (1/8")

1225-S [Short Threaded Bone Pins] **Pkg of 10**

Overall Length: 27mm

Collar-to-Tip Length: 13mm

Diameter: 3.2mm (1/8")

1025 [Sterilization Case]



MADE  
IN USA

FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.

**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

1.800.548.2362