


# Capsule Retractors

*Retractors for use in Bankart surgery*

The single prong retractor is commonly used when retracting on the inferior rim of the glenoid. The two and three-prong retractors are designed to be placed medially along the scapular neck to retract the anterior capsule and labrum.



#### PRODUCT NO'S:

T1008-01 [3 Prongs]  
Overall Length: 10" (25,4 cm)  
Prong Length: 30 mm

T1008 [2 Prongs]  
Overall Length: 10" (25,4 cm)  
Prong Length: 30 mm

T1009 [1 Prong]  
Overall Length: 10" (25,4 cm)  
Prong Length: 30 mm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

[www.innomed.net](http://www.innomed.net)

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**