


Cartilage Graspers

Helps to grasp and hold cartilage, tendons, soft tissues and loose bodies

Shaft allows for use in narrow spaces.



Saw Tooth Jaw




Saw Tooth Jaw with 6" (15,2 cm) Shaft.

PRODUCT NO:
1785 [Saw Teeth]
 Shaft Length: 6" (15,2 cm)
 Overall Length: 9.25" (23,5 cm)

MADE EXCLUSIVELY
 FOR INNOMED IN
 GERMANY



Shark Tooth Jaw



Shark Tooth Jaw with 8" (20,3 cm) Shaft.

Shark Tooth Jaw with 5" (12,7 cm) Shaft.

Shark tooth designed by Luis Ulloa,
 design modification by Michael Soudry, MD

PRODUCT NO'S:
1777 [5" with Shark Teeth]
 Shaft Length: 5" (12,7 cm)
 Overall Length: 8.25" (21 cm)
 Jaw Bite: 2 mm x 6.5 mm

1779 [8" with Shark Teeth]
 Shaft Length: 8" (20,3 cm)
 Overall Length: 11.25" (28,6 cm)

MADE EXCLUSIVELY
 FOR INNOMED IN
 GERMANY

www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
 info@innomed.net

912.236.0000 Phone
 912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
 info@innomed-europe.com

1.800.548.2362