

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2000 Innomed, Inc., All Rights Reserved.




Innomed
Instrument Division
103 Estus Drive, Savannah, GA 31404FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Cement Osteotome

The Cement Osteotome is designed to be inserted around the back of the tibia base to remove cement. The curve is congruent with most tibia bases. The osteotome is nitrate coated to help protect the implant surface.

Product No: 5220

Helps remove cement around the back of the tibia base



Cement Remover

The Cement Remover is designed with sharper faces to help remove unhardened cement around femoral and tibial knee components. The remover is nitrate coated to help protect implant surfaces.

Product No: 5230

Helps remove unhardened cement around femoral and tibial knee components