



Chandler Retractors

Used for retracting tissue away from the bone, and helpful for posterior exposure of the tibia in MIS surgery

Allows the surgeon to retract soft tissue away from the bone, and can be used for hip and knee surgery. The handle is contoured away from the field of view and working area. Available in three blade sizes: 5/8", 3/4" and 1".

The OrthoLucent™ version is made of a lightweight carbon fiber material, which is completely radiolucent, helps to prevent from scratching component surfaces, and can be steam sterilized.

PRODUCT NO.'S:	
3220-01 [5/8"] Overall Length: 9.125" Blade Width: 16mm	3220-04 [1"] Overall Length: 9.125" Blade Width: 25.4mm
3220-02 [3/4"] Overall Length: 9.125" Blade Width: 19mm	3220-02R [OrthoLucent™ 3/4"] Overall Length: 9.125" Blade Width: 19mm

OrthoLucent™ is a trademark of MedX Composites



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

INNOMED

1.800.548.2362