

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Chandler Retractor

Allows the surgeon to retract soft tissue away from the bone, and can be used for hip and knee surgery. The specially designed handle is contoured away from the surgeon's field of view and working area. Available in three blade sizes: 5/8" and 3/4".

Product No's: 3220-01 [5/8"]
3220-02 [3/4"]

Used for retracting tissue away from the bone

