


Charnley-type Self-Retaining Hip Retractor

Designed by
Perry Greene, MD

Product No: 1589
Overall Length: 12.25"
Length to Hinge: 9"
Leg Depth: 2.75"
Prong Dimensions: 12mm L x 4mm W x 3mm D
Opens to 7" Space between Legs



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2009 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

