

# Chen Diaphyseal Fracture Reduction Clamp

Designed by Franklin Chen, MD

*Designed to facilitate and maintain reduction of the internal fixation of diaphyseal and meta-diaphyseal fractures of long bones*

Works especially well with short oblique bones while providing room to implement the plate with this bone clamp still in place.



- ▶ Pivoting pads accommodate metaphyseal fractures
- ▶ The quick release enables adjustment without losing reduction
- ▶ Helps provide provisional reduction of diaphyseal fractures - humeral shaft fractures, tibial fractures


**PRODUCT NO:**

**1808**

Overall Length: 9.25" (23,5 cm)

Arm Downward Offset: 15 mm

Pad Dimensions: 1" x .375" (25,4 cm x 1 cm)



[www.innomed.net](http://www.innomed.net)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

**INNOMED**



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**