


Chung Weitlaner Retractor

Features long prongs allowing use of a small retractor in a small, but deep wound

Particularly useful for volar approach to the distal radius in large patients and for mini-open rotator cuff surgery.



Prong lengths of 25mm and 30mm available with either sharp or blunt tips

Designed by
Raymond Chung, MD

Product No's:
5065 BLUNT TIPS, 25MM BLADE DEPTH
5066 SHARP TIPS, 25MM BLADE DEPTH
5067 BLUNT TIPS, 30MM BLADE DEPTH
5068 SHARP TIPS, 30MM BLADE DEPTH

OVERALL LENGTH: 4.5"
OPENS TO MAXIMUM OF: 32MM



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

