


Cobra Retractors

A general purpose instrument for use around the femur and acetabulum

The OrthoLucent™ model is made of carbon fiber material—strong, lightweight, completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.



PRODUCT NO'S:	
6120 [Narrow] Overall Length: 11.75" Handle Length: 6.5" Blade Width: 16mm	6130-R [OrthoLucent™ Standard] Overall Length: 12" Handle Length: 7" Blade Width at Widest: 33mm
6120-04 [Extra Long Narrow] Overall Length: 15.5" Handle Length: 11" Blade Width: 16mm	6135 [Deep] Overall Length: 14.5" Blade Width at Widest: 33mm
6130 [Standard] Overall Length: 12" Handle Length: 7" Blade Width at Widest: 33mm	6140 [Wide] Overall Length: 11.75" Handle Length: 7" Blade Width at Widest: 56mm

OrthoLucent™ is a trademark of MedX Composites

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74



1.800.548.2362