

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

DMP Coring Instrument and Head Holder

A resected or allograft femoral head can be held with the Head Holder for reaming with the Coring Instrument. Cancellous bone is contained within the reamer head, which can be removed from the shaft for easy removal of the cancellous bone.

Designed by
DMP

Product No's:
8250 [Set]
8250-01 [Coring Instrument]
5250-02 [Reamer Head Replacement]
8260 [Head Holder]

Designed for handheld removal of cancellous bone from a resected femoral head

