

DMP Wire Tighteners

Used to tighten a cerclage wire around a bone

Can be used on wire up to 18 ga. diameter. T-Handle tightener is used to hand tighten a wire. Drill end tightener is used with a power drill or reamer—care must be taken to operate the power unit slowly to help eliminate the possibility of breaking the wires by over-tightening.

PRODUCT NO:

8729 [T-Handle]

Overall Length: 4.5"

Handle Width: 2.625"

End Diameter: 15mm

8630 [Drill Attachment End]

Overall Length: 4.625"

Designed by DMP



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED 

1.800.548.2362