



Exposure of the hip joint & removal of the femoral head



Acetabular exposure, reaming and cup insertion



Table Assembly/Elevator Hook for femoral access



Femoral broaching and stem insertion

Das/Seng Anterior Total Hip Instruments

Retractor set with included table-mounted controlled-release ratcheting elevator hook, specifically designed to help simplify anterior approach total hip arthroplasty

Designed by Amal Das, MD and Brian Seng, DO

Written Technique Available.



1 Posterior Femoral Neck Retractor

2 Anterior Femoral Neck Retractor

3 Anterior Acetabular Retractor

4 Posterior Acetabular Retractor

6 Femoral Calcar Retractor

7 Greater Trochanter Retractor

PRODUCT NO'S:

6221 [#1 - Posterior Femoral Neck Retractor]

1 Blade Width: 25mm
Blade Depth: 3"
Overall Length: 14"

6222 [#2 - Anterior Femoral Neck Retractor]

2 Blade Width: 31.5mm, 10mm @ Tip
Blade Depth: 4.5"
Overall Length: 15"

6223 [#3 - Anterior Acetabular Retractor]

3 Blade Width: 25mm
Blade Depth: 2.25"
Overall Length: 13.25"

6224 [#4 - Posterior Acetabular Retractor]

4 Blade Width: 25mm
Blade Depth: 2.75"
Overall Length: 14"

6226-EH [#5 - Elevator Hook]

5 Blade Width: 10mm
Blade Depth from T-Handle: 5.75"
Overall Length: 9.25"

6227 [#6 - Femoral Calcar Retractor]

6 Blade Width: 25mm
Blade Depth: 3.625"
Overall Length: 14"

6225 [#7 - Greater Trochanter Retractor]

7 Blade Width: 25mm
Blade Depth: 2.5"
Overall Length: 14.25"

6226-TA [#8 - Table Assembly]

8 This product number includes one 6226-EH Elevator Hook
Folds to approx: 21" x 5" x 5"



5 Elevator Hook

Table Assembly 8

MADE IN USA

FREE TRIAL ON MOST INSTRUMENTS

© 2012 Innomed, Inc.

INNOMEDI



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362