

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2002 Innomed, Inc., All Rights Reserved.

Innomed

Instrument Division
103 Estus Drive, Savannah, GA 31404




FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Modified Double Prong Acetabular Retractor with Center Prongs

Designed to retract the femur anteriorly during total hip arthroplasty. It is hooked over the anterior pelvic brim. Weights can be added to assist in exposure and to help hold the retractor in place.

Product No's: 6170
6175 [Narrow]



Retracts the femur anteriorly during total hip arthroplasty

