

Deep Curved Anterior Retractors

Deep retractors designed with a longer handle for use in Direct Anterior approach THR when added leverage and depth is needed for use in large patients


PRODUCT NO'S:	
7804-01D	[Deep Narrow] Overall Length: 15" (38,1 cm) Blade Depth: 7" (17,8 cm) Upper Blade Width: .6" (15 mm) Blade Width at Tip: .14" (3,5 mm)
7804-02D	[Deep Standard] Overall Length: 15" (38,1 cm) Blade Depth: 7" (17,8 cm) Upper Blade Width: .8" (21 mm) Blade Width at Tip: .14" (3,5 mm)
7804-03D	[Deep Wide] Overall Length: 15" (38,1 cm) Blade Depth: 7" (17,8 cm) Upper Blade Width: 1.44" (36 mm) Blade Width at Tip: .14" (3,5 mm)



Deep Narrow
#7804-01D

Deep Standard
#7804-02D

Deep Wide
#7804-03D



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.