


Spine/Trauma Deep Tissue Retractor

Designed to help maximize exposure with 90° arms and deep tissue blades

The retractor arms are available in configurations of 7 or 4 teeth.



PRODUCT NO'S:	
1862 [4 Teeth]	
Overall Length: 7.5" (19,1 cm)	
Handle-to-Bend Length: 6" (15,2 cm)	
Drop Depth: 3.25" (8,3 cm)	
Prongs: 1.5" Long x .75" Wide (38 mm x 19 mm)	
1863 [7 Teeth]	
Overall Length: 7.5" (19,1 cm)	
Handle-to-Bend Length: 6" (15,2 cm)	
Drop Depth: 3.25" (8,3 cm)	
Prongs: 1.5" Long x 1.375" Wide (38 mm x 35 mm)	

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

INNOMED

1.800.548.2362