


Desai Curette Osteotomes

Designed by Sarang Desai, DO

Designed to remove bone and cartilage, helpful for preparing joint surfaces for fusion, allowing easy removal of osteophytes and cartilage without having to switch instruments

The osteotome portion also can be used to "feather" the subchondral surface to expose bleeding bone. It is also useful in instances of obtaining autograft, as it can be used to create a bone window and then remove cancellous bone.



PRODUCT NO:

5242 [8 mm]
 Overall Length: 8.25" (21 cm)
 Handle Length: 4.25" (10,8 mm)
 Cup: 8 x 10 mm
 Osteotome Width: 6.5 mm
 Osteotome Length: 3 mm from edge of cup

Cup Size:

