




Digit Aligner

Radiolucent device designed to help reduce fractures of the small digits



Radiolucent



PRODUCT NO'S:

8 Hole Digit Aligner (Toe)

1190-8L [Large]

Outside Diameter: 1.840"

Inside Diameter: 1.515"

Overall Length: 1.240"

1190-8M [Medium]

Outside Diameter: 1.590"

Inside Diameter: 1.265"

Overall Length: 1.240"

1190-8S [Small]

Outside Diameter: 1.340"

Inside Diameter: 1.015"

Overall Length: 1.240"

12 Hole Digit Aligner (Finger)

1190-12L [Large]

Outside Diameter: 1.840"

Inside Diameter: 1.515"

Overall Length: 1.660"

1190-12M [Medium]

Outside Diameter: 1.590"

Inside Diameter: 1.265"

Overall Length: 1.660"

1190-12S [Small]

Outside Diameter: 1.340"

Inside Diameter: 1.015"

Overall Length: 1.660"

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362