

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2000 Innomed, Inc., All Rights Reserved.

**Innomed**  
Instrument Division  
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION


## Dozier Radiolucent Bennett Hip Fracture Retractor

The retractor is designed to be used in hip fractures with the advantage that the retractor can be kept in place while using image intensification or taking an x-ray. The handle is designed to be rotated to the right or left for surgeon preference. The retractor may be steam or gas sterilized.

*Designed by*  
John K. Dozier, MD

Product No: 6870

*Can be kept in place while using image intensification or taking an x-ray*



Shown with handle  
in both directions

