




# Duncan Metatarsal Clamp

Designed by Gregory S. Duncan, DPM

*Designed to clamp and hold an osteotomized metatarsal bone in the corrected position for fixation through the opening in the top of the clamp*



**PRODUCT NO:**

1638

Overall Length: 7"

Clamp Pads: 1.3" x .625"

Opening: 1" x .375"

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY

FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

1.800.548.2362