

Durham Bone Reduction Clamps

Designed by Alfred A. Durham, MD

Allows application of a bone plate without removing the reduction clamp

The large window directly above the jaws provide space to allow a bone plate to be slid into position without removing the clamp.

PRODUCT NO.'S:
3652 [Standard] Overall Length: 7" (17,8 cm)
3652-01 [Large with Speedlock] Overall Length: 9.25" (32,5 cm)



▶ The **standard clamp** is designed for medium size bones such as the fibula, ulna, and radius

Standard





▶ The **large clamp with speedlock** is designed for large bones such as the femur and tibia

Large with Speedlock



Large Window
Provides space to allow a bone plate to be slid into position without removing the clamp



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to Launch Our Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362