

Durkan Ratchet Bone Clamps


Designed by John Durkan, MD

Design of ratcheting mechanism allows for quick tightening and release around the bone

PRODUCT NO:

1867 [Large]
Overall Length: 8.625"
Jaw opens to: 3.5"

1868 [Small]
Overall Length: 8.5"
Jaw opens to: 3.75"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362