



Eickmann Knee Revision Set

Designed by Thomas Eickmann, MD

Used for total knee revision



PRODUCT NO'S:

5470-00 [Complete Set]

Individual Instruments:

5470-08 [8mm Chisel]

Osteotome Width: 8mm
Blade Length: 2.375"
Overall Length: 7.375"

5470-11 [11mm Chisel]

Osteotome Width: 11mm
Blade Length: 2.375"
Overall Length: 7.375"

5470-20 [20mm Chisel]

Osteotome Width: 20mm
Blade Length: 2.375"
Overall Length: 7.375"

5472-08 [8mm Offset Cement Removal Chisel]

Osteotome Dimensions: 8mm Wide x 12mm Long
Blade Length: 2.375"
Overall Length: 7.375"

5474-06 [6mm Notched Cement Removal Chisel]

Osteotome Width: 6mm
Blade Length: 2.625"
Overall Length: 7.375"

5475-08 [8mm Implant Remover]

Diameter: 8mm
Blade Length: 2.625"
Overall Length: 7.375"

5470-CASE [Case Only]

PROUDLY
MADE
IN THE
USA

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362