

Engel Unicondylar Minimally Invasive Cement Scrapers

Designed by Gerard A. Engel, MD

Designed for use in unicondylar minimally invasive knee surgery

Right and left design used to scrape cement from around and behind the implants.



PRODUCT NO'S:

4920-01 [Cement Scraper – Right]

Scraper Head: 5 mm x 9 mm

Overall Length: 8.5" (21,6 cm)

Handle Length: 3.5" (8,9 cm)

4920-02 [Cement Scraper – Left]

Scraper Head: 5 mm x 9 mm

Overall Length: 8.5" (21,6 cm)

Handle Length: 3.5" (8,9 cm)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

INNOMED