

# Evans Universal Carpal Tunnel Knife Guide

Designed by Peter J. Evans, MD, PhD

*Designed to protect the median nerve while providing a choice of grooved tracks for commercially available retrograde knives (that do not provide this feature) or for tenotomy scissors*

Allows for smooth advance of the blade or scissors to divide the transverse carpal ligament.  
Designed for a mini-open, non-endoscopic approach.

<b>PRODUCT NO:</b>
1128
Overall Length: 8
Blade Guide Widths: 2mm and 5mm



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362