


# Evans Modified Fukuda-type Retractors

Designed by Peter J. Evans, MD

*Designed to retract the humeral shaft posteriorly,  
helping to expose the glenoid surface*

Center groove allows the shaft of a reamer to fit more posteriorly.



#### PRODUCT NO'S:

##### 5180-N [Narrow]

Overall Length: 8.625" (21,9 cm)

Blade Width: 1" (2,54 cm)

Blade Depth: 3.75" (9,5 cm)

##### 5180-W [Wide]

Overall Length: 8.625" (21,9 cm)

Blade Width: 1.25" (31,7 mm)

Blade Depth: .75" (9,5 cm)

[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

# INNOMED

1.800.548.2362