

Evans Reverse Hohmann Retractor

Designed by Peter J. Evans, MD

*Designed for total shoulder arthroplasty
and open rotator cuff procedures*



Smaller size useful for retracting the deltoid muscle and other structures.

PRODUCT NO:

4547

Blade Width: Tapers from 30mm to 18mm

Blade Depth: 3"

Prong Width: 6mm

Overall Length: 8.5"

MADE
IN USA

FREE TRIAL ON MOST INSTRUMENTS

© 2012 Innomed, Inc.

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362