

# Extra Small Femur & Tibia Triangle

Designed by S.E. Fromm, MD  
& Kenneth Merriman, MD

*Used for femur and tibia positioning during nailing, repairs and fractures*

Supplied with an autoclavable silicone pad and velcro\* straps.  
The triangle is radiolucent and gas or steam sterilizable.



PRODUCT NO'S:	
2760-XS	Height: 8.5"
Replacement Parts:	
2760-P	[Silicone Pad]
8120-SP	[Straps] Package of 10

MADE IN USA

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

INNOMED

1.800.548.2362