

Femoral Canal Preparation Scrapers

Designed by Gary W. Bradley, MD

Used for femoral canal preparation in total hip surgery

Used to help remove excess bone in the inside of the femoral canal after broaching.
Replaces plastic brushes and eliminates potential plastic debris.


PRODUCT NO:

5460-02

Blade Width: 12mm

Overall Length: 17.25"

Handle Length: 4.25"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

INNOMED

1.800.548.2362